**Intended use**

Fully automated latex enhanced immunoassay for the detection of total immunoglobulin in human citrated plasma that reacts with Platelet factor 4 (PF4) when complexed to heparin or on the ACL TOP® Family. Heparin associated antibodies are commonly found in patients with heparin induced thrombocytopenia or thrombosis (HIT).

**Summary and principle**

Heparin-induced thrombocytopenia (HIT) is an immune complex-mediated disorder that can cause morbidity and mortality. In patients receiving heparin therapy, HIT is considered a primary disease that anticoagulant is administered to prevent thrombosis.<sup>1</sup> HIT is the major clinical event in HIT is an increased risk for venous and/or arterial thrombosis.<sup>2</sup> HIT is suspected when patients treated with unfractionated heparin (UFH) or low molecular weight heparin (LMWH) develop a thrombotic episode, which is often preceded by a period of days to weeks of normal platelet counts. Criteria of the American Society of Thrombosis and Hemostasis (ASTH) for HIT and delayed-onset HIT have also been reported.<sup>3</sup> Prompt diagnosis is very important since treatment must be suspended and alternative anticoagulants used in case of confirmed HIT. HIT, the immunological type, is caused by the development of platelet-activating antibodies, mostly directed against Platelet Factor 4 when complexed with Heparin (PF4-H).<sup>4</sup>

The current working hypothesis is that early exposure to PF4 complexes, perhaps induced by unfractionated heparin or anti-thrombin factors, leads to PF4-platelet and a T-cell-dependent sensitization. These antibodies are induced by heparin treatment.<sup>5</sup> Development of pathology is mainly associated with heparin-dependent antibodies of the IgG isotype.<sup>6</sup> However, the presence of only IgM and/or IgA was observed in some patients, including cases showing a thrombotic.

Calibrator and Control reagent after opening. Opened calibrator is stable 3 hours at 15-25°C on-board the ACL TOP Family.

**Instrument/test procedures**

Refer to the appropriate IL instrument's Operator's Manual and/or Application Manual for the complete assay procedure instructions.

**Specimen collection and preparation**

Nine parts of freshly drawn venous blood are collected into one part of trisodium citrate. Refer to CLSI Document H2-A5 for further instructions on specimen collection, handling and storage.<sup>7</sup>

Thaw frozen samples rapidly at 37°C and centrifuge before testing, if necessary. After thawing the assay must be performed within 7 hours.

**Additional reagents and control plasmas**

The following is not supplied with the kit and must be purchased separately.

**Hemosil HIT-Ab(PF4-H) Controls** Cat. No. 0020301200

**Factor Diluent** Cat. No. 0009757600

**Quality control**

Two levels of control are recommended for a complete quality control program.<sup>8</sup> Hemosil HIT-Ab(PF4-H) Control is a liquid dilution of monoclonal anti-PF4-Heparin in Tris buffer, containing bovine serum albumin, stabilizers and preservative.

**S** **Stabilizer** (Art. Nr. 0020301200) 2 vials x 0.8 mL of a solution of PIV-PVS complex (PF4 from human platelets complexed to PIV), in PBS buffer containing bovine serum albumin, stabilizers and preservative. Contains 0.02% Bronidox® as a preservative.

**C** **Calibrator** (Art. Nr. 0020301240) 2 flacons x 1 mL of monoclonal anti-PF4-Heparin in PBS buffer containing bovine serum albumin, stabilizers and preservative.

**Precautions and Warnings**

The Hemosil material used in this product was tested by FDA approved test methods and found non-reactive for Bovine Surface Antigen (BsAg), Anti-HCV and HIV 1/2 antibodies. Handle as if potentially infectious.

All reagents except the Complex contain less than 0.1% sodium azide that may form explosive azides in metal plumbing. Use proper disposal procedures.

**Latex Reagent, Stabilizer, Complex and Calibrator** contain bovine material. All donor animals were sourced to BSL-2 conditions. The cattle received annual and more mortal health inspection. All vitamins and preservatives were applied to reduce risk of infection and contagious material. However, the material should be treated as potentially infectious.

**Latex Reagent, Complex and Calibrator** are not classified as dangerous.

- Hazard class: none
- Risk phrases: none
- Safety phrases: none

**Stabilizer:**

- Hazard class: Xn (Harmful)
- Risk phrases: Harmful if swallowed (R22)
- Safety phrases: If swallowed, seek medical advice immediately and show this container or label (S46)

This product is For *in vitro* Diagnostic Use.

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