**HIT- Ab(PF4-H) - 0020301200**

**Intended use**

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**Summary and principle**

Heparin-induziert Thrombocytopenie (HIT) is an immune complex mediated disorder that can cause morbidity and mortality in patients receiving heparin therapy. HIT is considered a paradoxical disease in that anticoagulant therapy can be used to reverse thrombosis, however the major clinical event in HIT is an increased risk for venous and/or arterial thrombosis. HIT is suspected when patients treated with unfractionated heparin (UFH) or low molecular weight heparin (LMWH) show a decrease of platelets greater than 50% from baseline, typically between days 5 and 10 from the initiation of the anticoagulant treatment. Cases of early-onset HIT and delayed-onset HIT have also been reported. Prompt diagnosis is very important to ensure heparin treatment must be suspended and alternative anticoagulants used in cases of confirmed HIT. Type II HIT, the immunologic type, is caused by the development of platelet-activating antibodies, which are directed against heparin-anticoagulant complexes. These antibodies are the primary cause for inducing thrombosis both venous and arterial. The current working hypothesis is that early exposure to PF4 complexes, perhaps induced by endogenous non-heparin sources, leads to PF4 clustering and a T-cell-mediated antibody class switch of B-cells. Later in life, these B-cells become again transiently activated when PF4 clusters are induced by heparin treatment. Development of pathology is mainly associated with heparin-dependent antibodies of IgG type (IgG1). However, the presence of only IgM and/or IgA was observed in some patients, including cases showing a thrombotic complication.

PF4/H antibody testing combined with an appropriate clinical assessment has been proven to be useful as an aid in the management of HIT suspected patients. Particularly a negative result for a PF4/H antibody test can support the clinical decision to exclude the presence of HIT, and therefore continue heparin treatment. A weak positivity for PF4 antibodies may indicate that the antibodies are non-platelet activating, while a strong positivity may indicate a higher risk for HIT. In both cases, confirmation with a functional test is recommended. A clinical reassessment supported by laboratory data should be performed before confirmation or exclusion of the diagnosis. Hemosil HIT-AB<sub>PF4-H</sub> kit is a latex particle enhanced immunobimetric assay to detect total heparin associated antibodies found in HIT patients. A monoclonal antibody that mimics human HIT antibodies is coated on latex particles. In the presence of PF4 from human platelets complexed to polyvinyl sulfonate (PVS), and the patient sample, a competitive agglutination reaction occurs. The degree of agglutination is inversely proportional to the concentration of antibodies in the sample and is determined by measuring the decrease of transmitted light caused by the aggregates.

**Composition**

- [R]** **Latex Reagent** (Cat. No. 0020301210): 2 vials x 1.8 mL of a suspension of polystyrene latex particles coated with purified monoclonal anti-PF4-heparin in Tris buffer, containing bovine serum albumin, stabilizers and preservative.
- [S]** **Stabilizer** (Cat. No. 0020301220): 2 vials x 3.2 mL of PBS buffer containing bovine serum albumin, stabilizers and preservative.
- [Com]** **Complex** (Cat. No. 0020301230): 2 vials x 0.8 mL of a solution of PF4-PVS complex (PF4 from human platelets complexed to PVS) in PBS buffer containing bovine serum albumin, stabilizers and preservative.
- [C]** **Calibrator** (Cat. No. 0020301240): 2 vials x 1 mL of monoclonal anti-PF4-heparin in PBS buffer containing bovine serum albumin, stabilizers and preservative.

**Precautions and Warnings**

The human material used in this product was tested by FDA approved test methods found nonreactive for Hepatitis B Surface Antigen (HBsAg), Anti-HCV and HIV 1/2 antibodies. Handle as if potentially infectious. Do not mix reagents except the Complex contain sodium azide that may form explosive azides in metal plumbing. Use proper disposal procedures. Latex Reagent, Stabilizer, Complex and Calibrator contain bovine material. All donor animals were screened from BSE-free herds and are free from any transmissible spongiform encephalitis agent, and they were apparently free from infectious and contagious material. However, the material should be treated as potentially infectious.

- Hazard class: none
- Hazard phrases: none
- Precautionary phrases: none

**Preparation**

**Latex Reagent:** Gently invert several times to mix before use. Do not shake. Avoid foam formation. **Stabilizer:** Gently mix several times before use. Do not shake. Avoid foam formation. **Complex:** Gently mix several times before use. Do not shake. Avoid foam formation. **Calibrator:** Gently mix several times before use. Do not shake. Avoid foam formation. Note: Bubbles on top of the liquids may interfere with the instrument's liquid sensors.

**ENGLISH - Insert revision 11/2014**

**Limitations/interfering substances**

HIT-AB<sub>PF4-H</sub> results on the ACL TOP Family are not affected by hemoglobin up to 495 mg/dL, bilirubin up to 18 mg/dL, triglycerides (Intralipid®) up to 375 mg/dL, rheumatoid factor up to 10 IU/ml, and human anti-mouse antibodies (HAMA) up to 1 µg/ml. Samples from non-HIT suspected, heparin treated patients (either UFH or LMWH) were tested with the Hemosil HIT-AB<sub>PF4-H</sub>. For LMWH, all twenty results were below the recommended cut-off in a sample population with a heparin concentration range between 0.25-1.08 IU/ml. For UFH, this cut-off of twenty-two were slightly above the cut-off in a sample population with a heparin concentration range between 0.15-74 IU/ml. However, there was no dose-response correlation between heparin concentration and UFH concentration.

**Expected values**

The presence of PF4/H antibodies in a normal population is not expected. A set of studies was performed to assess the optimal cut-off of Hemosil HIT-AB<sub>PF4-H</sub> on the ACL TOP Family systems as described below. 151 plasma samples from apparently healthy individuals (blood bank donors) were tested. The 95% reference interval was 0.6-10 U/ml. In addition, a population of 51 samples from heparin-exposed patients (UFH and LMWH) with no clinical suspicion of HIT was tested. The 95% reference interval was 0-12 U/ml.

Finally, a method comparison with the Serotonin Release Assay (SRA) on 66 plasma samples from HIT suspected patients (31 SRA positive and 35 SRA negative) indicated that the optimal cut-off determined by ROC analysis was 1.0 U/ml (92.4% Agreement, CI% = 83.2% - 97.5%). Based on these studies it has been determined that on heparin treated patient samples, Hemosil HIT-AB<sub>PF4-H</sub> results equal or higher than 1.0 U/ml may indicate the presence of HIT antibodies. The assay results should be used with other information, including the clinical context, in forming a diagnosis.

**Performance characteristics**

**Precision:** Within run (repeatability) and total (within device) precision was assessed over multiple runs. The Hemosil HIT-AB<sub>PF4-H</sub> assay results are summarized below.

ACL TOP Family	Mean (U/ml)	CV% (Within run)	CV% (Total run)
Low HIT-AB <sub>PF4-H</sub> Control	0.7	0.09 (SD)	0.11 (SD)
High HIT-AB <sub>PF4-H</sub> Control	3.5	4.9	6.1
Weakly Positive HIT-AB <sub>PF4-H</sub> sample	1.6	4.9	8.1
High HIT-AB <sub>PF4-H</sub> plasma sample	5.2	2.8	3.5
Very High HIT-AB <sub>PF4-H</sub> plasma sample	19.1	2.8	3.5

**Method Comparison**

An in-house method comparison was performed on 257 frozen samples from HIT suspected patients using a commercially available ELISA test. The results summarized below are based on a cut-off of 1.0 U/ml for the Hemosil HIT-AB<sub>PF4-H</sub> assay.

	ELISA Results	+	-	14	148
<b>N</b>	<b>% Co-Positivity (95% CI)</b>	<b>% Co-Negativity (95% CI)</b>	<b>% Agreement (95% CI)</b>		
257	85.9% (77.4% - 92.0%)	93.7% (88.7% - 96.9%)	90.7% (86.4% - 93.9%)		

	ELISA Results	+	-	14	148
<b>N</b>	<b>% Co-Positivity (95% CI)</b>	<b>% Co-Negativity (95% CI)</b>	<b>% Agreement (95% CI)</b>		
257	85.9% (77.4% - 92.0%)	93.7% (88.7% - 96.9%)	90.7% (86.4% - 93.9%)		

**Clinical performance**

A multicenter study was performed at three hospitals on 414 samples comparing with a commercially available ELISA assay. The samples were from patients exposed to heparin who showed HIT related symptoms and hence were suspected to suffer from HIT.

	ELISA Results	+	-	33	313
<b>Predictive Device</b>	<b>% Co-Positivity (95% CI)</b>	<b>% Co-Negativity (95% CI)</b>	<b>% Agreement (95% CI)</b>		
ELISA assay	414	60.2% (48.9% - 70.8%)	94.6% (91.5% - 96.7%)	87.7% (84.1% - 90.7%)	

The precision, method comparison and clinical study results were obtained using specific lots of reagents and controls.

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