

The Fritsma Factor

YOUR INTERACTIVE HEMOSTASIS RESOURCE



Improving Acute Care Using Coagulation Mixing Studies

George A Fritsma MS, MLS

The Fritsma Factor,

Your interactive Hemostasis ResourceSM

Sponsored by Precision BioLogic

Dartmouth, Nova Scotia

george@fritsmafactor.com

www.fritsmafactor.com



Acute Care Mixing Studies Bottom Line at the Start (BLAST)

The participant...

- Lists the clinical applications for PT and PTT mixing studies.
- Lists the steps to perform a PTT mixing study.
- Reviews the technical demands of mixing studies.
- Explains why the mixing study is an acute care assay.
- Correlates mixing study results with lupus anticoagulant and specific inhibitor testing.

Mixing Study: An Acute Care Assay

Differentiates a coagulation factor deficiency from a specific inhibitor or from a lupus anticoagulant.



Kershaw GK, Orellana D. Mixing tests: diagnostic aides in the investigation of prolonged prothrombin times and activated partial thromboplastin times. *Semin Thrombos Hemost* 2013;39:283–90.

32-year-old Female Mucocutaneous Bleeding

- Two weeks post-partum
- Easy bruising, frequent nosebleeds, vaginal bleeding



Mucocutaneous Bleeding 32-yo Female, 2 Weeks Post-partum

Assay	Patient	RI
HGB	11.8 g/dL	12–15 g/dL
PT	12.4 s	9.8–12.6 s
PTT (APTT)	42.5 s	25–35 s
PLT count	310,000/ μ L	250–450,000/ μ L
Fibrinogen	320 mg/dL	220–498 mg/dL
Isolated, prolonged PTT—response? 1:1 PTT mix		

Isolated Prolonged PTT: Differential

- Could be nothing: 5% of normals exceed limit ($\pm 2SD$)
- Outpatient: DOAC, especially DTI: dabigatran?
- Inpatient: unreported UFH
- Congenital single factor deficiency: VIII, IX, or XI, hemophilia A, B, or C with bleeding; VWD
- Congenital FXII, PK, or HMWK without bleeding
- Acquired FVIII or IX inhibitor with severe bleeding
 - Auto-anti-factor VIII, “Acquired hemophilia”
- Lupus anticoagulant (LA, LAC)

DOACs: "The Girls"



Trixi

Fonda



Riva



Pixi



Edi



Dabi



Direct Oral Anticoagulants (DOACs)

- Oral direct anti-Xa anticoagulants: prolong the PT
 - Rivaroxaban; Johnson and Johnson, Xarelto
 - Apixaban: Bristol-Meyers Squibb, Eliquis
 - Edoxaban; Daiichi-Sankyo, Savaysa
 - Betrixaban; The Medicines Company, Bevyxxa
- Oral direct thrombin inhibitor (DTI): prolongs the PTT
 - Dabigatran; Boehringer-Engelheim, Pradaxa
- Intravenous DTI
 - Argatroban; ACOVA, generic available, PTT
 - Bivalirudin; Angiomax, PTT



Mostly




Mostly

PTT Mixing Study: Cheap and Simple


- Mix patient plasma 1:1 with pooled normal plasma (NP) and perform immediate PTT on mixture
- No correction: If PTT of 1:1 mix is $>10\%$ longer than NP PTT
 - Lupus anticoagulant (LAC)
- Correction: If PTT of 1:1 mix is $\leq 10\%$ longer than NP PTT
 - Factor deficiency?
 - Specific inhibitor (anti-FVIII)?, requires 37°C incubation

PTT Mixing Study



The Fritsma Factor

YOUR INTERACTIVE HEMOSTASIS RESOURCE




PTT Mixing Study Using 10% Limit

PTT
≤33 s: Correction
>33 s: No correction

PTT Mix: Why Does This Work?

- Hypothetical 20% F VIII level prolongs PTT
 - PTT rgts are calibrated to prolong at 30–40% FVIII, IX, XI
- Add NP with established 100% factor level
 - 1:1 mix, average of 100% and 20% = 60%, PTT corrects
- Hypothetical anti-FVIII or lupus anticoagulant
 - With typical avidity, retains its ability to prolong the mix



PT and PTT Results can Predict Inherited Single-factor Coagulopathies

PT	PTT	Single Factor Deficiency
Long	Normal	VII
Long	Long	X, V, II, and fibrinogen ¹
Normal	Long	VIII, IX, XI ²

¹PT & PTT prolonged only when fibrinogen is <100 mg/dL, perform fibrinogen assay (fibrinogen assay often added to initial screen)

²Contact factor deficiencies XII (1–3% prevalence), prekallikrein (PK, Fletcher), or high molecular weight kininogen (HMWK, Fitzgerald) also prolong PTT results, but not associated with bleeding

The Fritsma Factor


YOUR INTERACTIVE HEMOSTASIS RESOURCE

PTT reagent: Ca^{++} , particulate negative charge activator, phosphatidyl serine; test prolonged by XII, PK, HMWK, XI, IX, VIII, X, V, prothrombin, Fg deficiency, UFH, DTI, LAC

**Extrinsic
(Tissue Factor)**




PT reagent: tissue factor, Ca^{++} , phosphatidyl serine; prolonged by VII, X, V, prothrombin, Fg deficiency; coumadin Rx, **direct anti-Xa Rx**, LAC?



Intrinsic (Contact)

Figure courtesy of Margaret G. Fritsma, Rodak's Hematology, 5th Edition, 2015

Common




1:1 PTT Mix with Incubation

- PTT of immediate mix $\leq 10\%$ longer than NP
 - Correction: factor deficiency? But first...
 - Incubate 1:1 mix 1–2 hours and repeat
- Correction after incubated mix = factor deficiency
- No post-inc correction: PTT remains $>10\%$ above NP
 - Specific inhibitor such as anti-FVIII;
 - IgG₄: Temp dependent, requires incubation
- NP must be incubated. Why?

Thom J, Ivey L, Eikelboom J. Normal plasma mixing studies in the laboratory diagnosis of lupus anticoagulant. *J Thromb Haemost* 2003;1:2689–91

37°C Incubated 1:1 PTT Mix




Incubated Mixing Study Result 32-yo Female, 2 Weeks Post-partum

Assay	Result		Comment
PTT	42.5 s	RI: 25–35 s	Confirms previous PTT
PTT/NP1:1 mix immediate	32.1 s	Limit: 33 s	NP: 30s
PTT/NP1:1 mix 1–2 h at 37°C	37.3 s	Incubated NP: 35s, limit 38.5 s (10%)	

Conclusion: both immediate *and* incubated mix PTTs correct, suspect factor deficiency, arrange for factor assays and von Willebrand disease profile

Factor Assay Results

32-yo Female, 2 Weeks Post-partum

Factor	Result	RI	Comment
VIII	32%	50–150%	 ?
IX	92%		
XI	131%		
XII	113%		
HMWK	ND	65–135%	XII, HMWK & PK deficiency not associated with bleeding
PK			

Expanded PT and PTT Results

Deficiency/Disorder	PT	PTT	PT Mix	PTT Mix
VWF ¹ , VIII, IX, XI	N	Long	ND	Correc
VII	Long	N	Correc	ND
FG, II, V, X	Long	Long	Correc	Correc
Liver Dis, VK Def, DIC, Amyloid ²				
Nephrotic Syndrome	N	Long	ND	Correc
LAC, VIII, IX Inh	N	Long	ND	No Cor
DOACs ³	Long	Long	No Cor	No Cor
UFH	N	Long	ND	No Cor

1: if VIII < 35 U/dL; 2: if X deficient; 3: Anti-Xa prolongs PT, DTI prolongs PTT

52-yo Athletic Female

Screen Prior to Hip Replacement Surgery

Test	Result	RI
HGB	14.1 g/dL	12–15 g/dL
PT	11.2 s	9.8–12.6 s
PTT	58 s	25–35 s
PLT	170,000/ μ L	150–400,000/ μ L
Fibrinogen	410 mg/dL	220–498 mg/dL


She reports no bleeding or bruising, no thrombosis

Again: Isolated Prolonged PTT: Differential

- Could be nothing: 5% of normals exceed limit ($\pm 2SD$)
- Outpatient: DOAC, esp DTI
- Inpatient: unreported UFH
- Congenital single factor deficiency: VIII, IX, or XI, hemophilia A, B, or C with bleeding; VWD
- Congenital FXII, PK, or HMWK without bleeding
- Acquired FVIII or IX inhibitor with severe bleeding
 - Auto-anti-factor VIII, “Acquired hemophilia”
- Lupus anticoagulant (LA, LAC)

52-yo Female PTT Mixing Study

Test	Result	Comment
TT	17 s	RI: < 21 s, rules out dabigatran
PTT	58 s	RI: 25–35 s
PTT NP	28 s	Correction if < 30.8 s (10%)
1:1 mix	35 s	25% longer than NP = no correction
What is the next step?		



PTT Reagent LAC Response

Fritsma Factor 2011–12 Survey: LAC Sensitivity

LAC Response	FF	CAP 2011	Examples
High	10%	46%	Siemens Actin FSL HemosIL aPTT-SP
Intermediate	43%	47%	Stago STA-PTT-LA HemosIL SynthASil
Low	18%	3%	Siemens Actin FS Stago CK–Prest
Don't know	29%	Other: 4%	


Fritsma GA, Dembitzer FR, Randhawa A, et al. Recommendations for appropriate activated partial thromboplastin time reagent selection and utilization. *Am J Clin Pathol* 2012; 137: 904–8.

Variation in PTT Reagent Response to LAC

LAC Pos Pat	Triniclot APTT S	Triniclot APTT HS	Actin FSL	Actin FS	PTT-LA
1	34	61	65	31	132
2	62	74	74	34	138
3	88	186	164	47	>300
4	234	64	69	30	253
5	60	171	136	66	258
Limit	154	35	40	36	54

Kershaw GK, Orellana D. Mixing tests: diagnostic aides in the investigation of prolonged prothrombin times and activated partial thromboplastin times. Semin Thrombos Hemost 2013;39:283–90.

Acute Care PTT Mixing Study Algorithm





The Fritsma Factor

YOUR INTERACTIVE HEMOSTASIS RESOURCE

Mixing Study Considerations

Specimen Management

Preameritics

Specimen Collection Factors

- Hemolysis: reject
- Prolonged tourniquet application raises VWF concentration
- Lipemia and icterus: use electro-mechanical coagulometer
- Wrong anticoagulant: EDTA (lavender), heparin
 - Run a potassium and a lithium
- Underfilled tube: AC proportionally increased
- Combining two tubes: AC proportionally increased
- Partially activated or clotted specimen: reject
- Specimen must be platelet-poor, $<10,000/\mu\text{L}$
 - Platelets release platelet factor 4, neutralizes heparin, also factor V

More Preanalytical Variables

- Start within 2 h to avoid in vitro prolongation
 - FV and FVIII deteriorate
- Coumadin, anti-Xa rivaroxaban, apixaban, edoxaban, betrixaban prolong PTT, PT
- Dabigatran, IV DTIs, and UFH prolong PT, PTT Mostly
- Heparinase neutralizes ≤ 1 unit/mL UFH
- 15% of anti-FVIII inhibitors detected in immediate mix
- 15% of LAs require incubation
- Weak LAs may be missed in 1:1 mix
 - Select a more LA-responsive PTT reagent or prepare a 4:1 mix

Normal Plasma Source?

- Home brew: pool ~20 normal plasmas, male \cong female
 - Ensure plasma is platelet-poor; $< 10,000/\mu\text{L}$; PTT \cong mean of RI
 - Ensure NP has ~100% of all factors, especially VIII, IX, and XI
 - Elevated FVIII causes false negative results (patient specimens)
 - Screen each for LA, specific factor inhibitors. HBV, HCV, HIV
 - Aliquot and freeze
- Or purchase commercial plasma
 - GMP & frozen meets all criteria: Precision BioLogic, George King Biomedical
 - Lyophilized plasma acceptable when validated
 - Processed with stabilizers

Clinical and Laboratory Standards Institute. One-stage prothrombin time (PT) test and activated partial thromboplastin time test (APTT) approved guideline—second edition. CLSI Document H47-A2. CLSI, Wayne PA. 2008.

What Limit Defines Correction?

No Consensus; Fritsma Factor 2015 Quick Question Answers

- Limits based on a fixed PTT value such as reference interval
 - 1:1 mix within RI upper limit (95% or 99% confidence interval, 39%)
 - 1:1 mix within RI upper limit + 5 seconds (8%)
- Limits based on the pooled normal plasma PTT value
 - 1:1 mix within NP PTT value + 5 seconds (14%)
 - 1:1 mix within NP PTT + 10% (32%)
- Rosner or Chang limit formula using patient, NP, and 1:1 mix results
 - Rosner formula produces a ratio
 - Chang's formula produces % deviation, requires incubation of patient plasma
- Other (7%): some combination of RI and Rosner, or...
 - A dedicated RI using ~20 healthy volunteer samples tested in 1:1 mixes. then computing 95th or 99th percentile.

Chang Index: Limit Based on % Correction

$$\% \text{ Correction} = \frac{\text{Patient PTT} - 1:1 \text{ mix PTT}}{\text{Patient PTT} - \text{NP PTT}} \times 100$$

$$\% \text{ Correction} = \frac{42.5 - 32.1 = 10.4}{42.5 - 30 = 12.5} = 0.83 = 83\%$$

Factor Deficiency = $\geq 75\%$

Inhibitor = $< 75\%$

Incubate patient plasma, NP, and mix

% Correction verified by local laboratory

Chang SH, Tillema V, Scherr D. A "percent correction" formula for evaluation of mixing studies. Am J Clin Pathol 2002;117:62-73.

Rosner Index

$$\text{Rosner Index} = \frac{1:1 \text{ mix PTT} - \text{NP PTT}}{\text{Patient PTT}} \times 100$$

$$\text{Rosner Index} = \frac{32.1 - 30}{42.5} \times 100 = 4.9$$

Inhibitor

≥ 11

Have to incubate patient specimen

Correction

< 11

Slightly more conservative than 10

Rosner index validated by local laboratory

Rosner E, Pazner R, Lusky A, Modan M, Many A. Detection and quantitative evaluation of lupus circulating anticoagulant activity. *Thromb Haemost* 1987; 57: 144-147.

Limit Levels and Correct Interpretation

Method	Result	Factor Def	Inhibitor
Mix secs – NP secs	<4 s	88%	100%
	<8 s	100%	77%
<u>Mix secs</u> NP secs	<1.1	81%	100%
	<1.2	98%	82%
Rosner Index	<11	100%	82%
	<15	100%	61%
Chang Index	>72%	100%	59%

The Fritsma Factor

YOUR INTERACTIVE HEMOSTASIS RESOURCE

59-yo Male Former Hockey Player

Total knee replacement preop labs




59-yo Male Former Hockey Player

Screen Prior to Knee Replacement Surgery

Test	Result	RI
HGB	14.8 g/dL	12–15 g/dL
PT	11.2 s	9.8–12.6 s
PTT	38 s	25–35 s
PLT	310,000/ μ L	150–400,000/ μ L
Fibrinogen	390 mg/dL	220–498 mg/dL
Patient reports no bleeding or bruising, no thrombosis		

When to Perform Mixing Study?

- Any PTT $>$ RI upper limit
- Any PTT $>$ RI upper limit + 5 seconds
- Any PTT $>$ RI upper limit with consult
- Is patient bleeding or clotting? 
- Anticoagulant?
- Possible mild coagulopathy?
- Possible “weak” LA?

59-yo Male

Former Hockey Player: TKR

Test	Result	Comment
TT	17 s	RI: < 21 s, rules out dabigatran
Chromo anti-Xa	0.0	0.4–0.7 U/mL, r/o direct anti-Xa
PTT	38 s	RI: 25–35 s
NP	31 s	Correction if < 34.1 s (10%)
1:1 mix	35 s	Correction? No correction?
What is the next step?		

59-yo Male Former Hockey Player Clinical Consult


- Consult: if no medical conditions, go on to surgery
- If he has had prior VTE: DVT or PE...
 - Perform mix using 4:1 patient plasma to NP
 - Or choose PTT reagent that is LA-sensitive
- If anatomic bleeding, test for inherited VWD or FVIII, FIX, FXI deficiency
- Or acquired vitamin K deficiency, liver disease, malignancy, renal insufficiency

The Fritsma Factor

YOUR INTERACTIVE HEMOSTASIS RESOURCE

2 YO Hemophilic Boy

Bleed into knee and ankle



2-yo Hemophiliac


Test	Result	RI
HGB	11.8 g/dL	9.6–15.6 g/dL
PT	11.2 s	9.8–12.6 s
PTT	65 s	25–35 s
PLT	310,000/ μ L	150–400,000/ μ L
Fibrinogen	390 mg/dL	220–498 mg/dL
Inflamed, swollen knee and ankle		

Mixing Study Result 2-yo Hemophilic Boy

Assay	Result	RI/NP	Comment
PTT	65 s	25–35 s	Confirms previous PTT
PTT/NP 1:1 mix immediate	33.5 s	NP 30 s	Correction (ambiguous)
PTT/NP 1:1 mix 2 h at 37°C	47.9 s	NP 35 s	NP is incubated alone and with mix
Conclusion: Anti-FVIII inhibitor			

Factor VIII Assay

- Dilute plasma 1:10
- Add factor VIII-depleted reagent plasma 1:1
- Add PTT reagent, incubate 3 minutes
- Add CaCl_2 , record interval to clot formation
- Compare result in seconds to calibration curve



Factor VIII Assay Dilutions Parallelism Indicates No Inhibitor

Plasma Dilution	Seconds	Raw Factor VIII Activity	Computed Factor VIII Activity (\times dilution)
1:10 "undil"	90 s	20%	20%
1:20	104 s	10%	20% (parallel)*
1:40	107 s	5%	20% (parallel)
1:80	110 s	2.5%	20% (parallel)

* $<10\%$ difference from undiluted indicates parallelism, no inhibitor

FVIII Assay Dilutions non-Parallelism Indicates Inhibitor

Plasma Dilution	Seconds	Raw Factor VIII Activity	Computed Factor VIII Activity (\times dilution)*
1:10 "undil"	80 s	10%	10%
1:20	93 s	8%	16%
1:40	107 s	5%	20%
1:80	108 s	4%	32%

* $>10\%$ difference from undiluted, rising = non-parallel, implies inhibitor

Kasper CK. Laboratory diagnosis of factor VIII inhibitors. In Kessler C, Garvey MB, Green D, Kasper C, Lusher J. Acquired Hemophilia 2nd Edition. Excerpta Medica 1995

The Fritsma Factor

YOUR INTERACTIVE HEMOSTASIS RESOURCE

55 YO Male with Atrial Fibrillation



55-yo Male with Atrial Fibrillation

Test	Result	RI
HGB	13.8 g/dL	12–15 g/dL
PT	17.2 s	9.8–12.6 s
PTT	159 s	25–35 s
PLT	310,000/ μ L	150–400,000/ μ L
Fibrinogen	20 mg/dL	220–498 mg/dL

55-yo Male with Afib Mixing Study

Assay	Result	RI
PTT	159 s	25–35 s
TT	> 150 s	< 21 s
PTT/NP 1:1 mix immediate	78 s	NP 30 s
PT/NP 1:1 mix immediate	15.2 s	NP 12 s
What do you recommend?		

If the PT is Prolonged

- Congenital deficiencies of II, V, VII, or X
 - PT and PTT long: II, V, X
 - PT only: VII, skip mixing and go to factor assay
 - Prevalence: 500,000–1:2,000,000
- Vitamin K deficiency: des-carboxy II, VII, and X
- Liver disease: PT prolongs before PTT due to des-carboxy II, VII, and X, reduced factor V

Develop Mixing Study Reliability

- PTT reagent sensitivities to factors and to LA
 - Activator: ellagic acid, silica, kaolin, celite
 - 30–40% FVIII, FIX, FXI
 - Intermediate sensitivity to LA
- NP consistency: ~100% activity for all factors
- Consultation for equivocal patient results
- Employ consistent correction limit

DIY Local Mixing Studies—Why?

- Unexpected isolated prolonged PTT or PT may require immediate therapy
- Local results may immediately direct therapy
- Delayed specimen may deteriorate
- Forward mixing study results to ref lab to direct follow-up, for instance, LAC profile or Bethesda titer



The Fritsma Factor

YOUR INTERACTIVE HEMOSTASIS RESOURCE

Lupus Anticoagulant

Advanced Mixing Study


Lupus Anticoagulant (LA)

- Requires two initial assays
 - Reduced phospholipid concentration PTT (PTT-LA[®])
 - Dilute Russell viper venom time (DRVVT, DVVtest[®])
- LAC present if either PTT or DRVVT is...
 - Prolonged with no correction upon mixing with NP
 - Corrected by high phospholipid reagent
 - PTT neutralization by increased phospholipid (Sta-Clot LA[®])
 - DRVVT confirmation by phospholipid reagent (DVVconfirm[®])
- Other causes ruled out...
 - Acquired multiple factor deficiencies
 - Factor-specific inhibitors

Based on **GOBSAT** research: Pengo V, Tripodi A, Reber G, et al. Update of the guidelines of lupus anticoagulant detection. J Thromb Haemost 2009;7:1737–40

Coag Mechanism (Again)

**Extrinsic
(Tissue Factor)**



Intrinsic (Contact)

DRVVT reagent:
Russell viper venom,
Ca⁺⁺, and
phospholipid activate
X, prolonged by
deficiency of X, V, II,
or fibrinogen;
coumadin Rx, LA

Common





Figure courtesy of
Margaret G. Fritsma,
Rodak's Hematology,
5th Edition, 2015

Lupus Anticoagulant-Sensitive Partial Thromboplastin Time (PTT-LA[®]) and confirmatory StaClot LA[®]




1. If the DRVVT is normal, proceed to the Hepzyme[®] step. If the DRVVT is prolonged, Hepzyme[®] is unnecessary.
2. If the StaClot LA[®] is positive and the DRVVT is negative, assay FVIII to determine if the positive StaClot LA[®] is due to a FVIII inhibitor.

The Fritsma Factor

YOUR INTERACTIVE HEMOSTASIS RESOURCE

Dilute Russell Viper Venom Time (DRVVT) and Confirmation
DRVVT First step = LA-Check; Second step = LA-Confirm



LAC By The Numbers

PTT-LA RI: 34–50 s

When PTT-LA is >50 s R/O heparin using thrombin time

Thrombin time RI: 17–21 s

When Staclot-LA Δ is ≥ 8 seconds, LAC is confirmed

DRVVT: LA Check RI : 30.9–41.5 s

LA-Sure ratio reference interval: 0.87–1.19

When LA-Sure ratio is ≥ 1.2 , LAC is confirmed


If Staclot-LA confirms LAC but DRVVT does not, assay factor VIII
If FVIII $< 30\%$, use Nijmegen Bethesda assay for FVIII inhibitor

The “LAC Cofactor” Effect

- Initial PTT = 48 s, RI 25–35; 1:1 mix *prolongs* to 54 s
- LAC binds IIa, slows clot formation, NP in mix adds IIa?
 - Or placental annexin V?
 - Do we really know?
- Magrath M. Lupus cofactor phenomenon. Letter J Clin Pathol 1990;42:264.
- Rand JH, Wu XX, Andree HA, et al. Antiphospholipid antibodies accelerate plasma coagulation by inhibiting annexin-V binding to phospholipids: a "lupus procoagulant" phenomenon. Blood. 1998;92:1652–60.
- Clyne LP. Plasma requirement for expression of lupus-like anticoagulant. Folia Haematologica int Ma Klin Morphol Blutforsch 1986;113:841

Anti-phospholipid Antibodies (APLs)

Antibodies that bind protein-anionic phospholipid complexes



Immunoassays for low-incidence autoantibodies

- Anti-phosphatidylserine
- Anti-phosphatidylcholine
- Anti-phosphatidylethanolamine
- Anti-prothrombin
- Anti-annexin

Immunoassays for anti-cardiolipin antibody

Clot-based assays for lupus anticoagulant

Immunoassays for anti-β₂-glycoprotein 1 antibody


Nigel Harris and Aziz Gharavi

APLs: Antibodies to Protein-Anionic Phospholipid Complexes

- Lupus anticoagulant (LA)
 - Defined by clot-based assays previously described
- Anti-cardiolipin antibody (ACL)
 - Defined by cardiolipin-based immunoassay
- Anionic phospholipid-binding proteins
 - β_2 glycoprotein 1 (target for standard immunoassay)
 - Annexin V (placental)
 - Prothrombin

} Minor targets

β_2 GP1 54KD single-chain glycosylated protein with five *sushi* domains binds anionic phospholipids



Anti- β_2 GP1 Vs. ACL

- ACL inter-laboratory standardization weak
- Pure β_2 glycoprotein 1 target antigen
 - No phospholipids in assay system
 - Antibodies associated with thrombotic events
- Ten β_2 GP1 kits
 - Pos agreement among kits IgG: 12/22, IgM: 5/22
 - No standard units nor cut-off values


Reber G, Tincani A, Sanmarco M, et al. Proposals for the measurement of anti- β_2 -glycoprotein I antibodies. Standardization group of the European Forum on Antiphospholipid Antibodies. *J Thromb Haemost* 2004;2:1860–2.

Enzyme Immunoassay Technology

- Anti- β_2 GP1 or ACL
- Report ACL as anticardiolipin antibody IgG or IgM units (APU)
- Report β_2 GP1 as β_2 GPI IgG or IgM units

Schematic Representation of Solid Phase ELISA for Detection of Anti-Beta 2 Glycoprotein 1 Antibodies

Antigen coated
microtiter wells




Step 1: Incubate 15 minutes with diluted patient, control, or calibrator sera.



Wash four times (PBS)

Step 2: Incubate 15 minutes with HRP-conjugated anti-human IgG, IgM, or IgA.



Wash four times (PBS)

Step 3: Incubate 10 minutes with chromogenic substrate (TMB); stop reaction; measure absorbance.



Spectrophotometer reading at 450 nm

Alloimmune APLs

- Transient, no thrombosis
- 2–4% of unselected individuals, rises with age
- Associated with infections, cancer, syphilis, HIV, drugs
 - Bacterial, viral, fungal, and parasitic infections
- Immunoassay: ACL or β_2 GPI positive
- Require confirmation after 12 weeks
- Clot-based assays variable

Diott JS, Belter M, Ercolano E, et al. Evaluating the role of a new dilute prothrombin time assay ACTICLOT dPT in lupus anticoagulant testing. *Blood* 2007;110:1061A

Autoimmune APLs

- SLE, RA, AITP, AIHA, other autoimmune disorders
- Primary anti-phospholipid syndrome
- Persistent LAC, APL positive
- Repeated spontaneous abortion related to ACL
- 30% associated with thrombosis
- Reported as chronic APL if repeat positive after 12 weeks

Anti-phospholipid Syndrome (APS)

- Arterial thrombosis
 - Transient ischemic attacks (TIAs), strokes
 - Acute myocardial infarctions (AMIs)
- Venous thromboembolic disease (VTE)
 - Deep venous thrombosis (DVT)
 - Pulmonary emboli (PE)
- Recurrent 2nd trimester fetal loss
 - Infertility, prematurity, intrauterine growth retardation, low birth weight
- Occasionally thrombocytopenia

Clinical Conditions Associated with APL

- Autoimmune diseases: SLE, RA, Sjogren syndrome, ITP, autoimmune hemolytic anemia
- Neoplasms: Hodgkin disease, myelofibrosis, epithelial malignancies
- Lymphoproliferative disorders: hairy cell leukemia, malignant lymphoma, Waldenstrom macroglobulinemia
- Drugs: chlorpromazine, procainamide, hydralazine, quinidine, antibiotics—penicillins, cephalosporins, phenytoin
- Infections: bacterial, protozoan (*P carinii*), viral

APL Consequences

- APLs do not increase risk of surgical bleeding.
- APL increases risk of post-surgical thrombosis.
- The PTT cannot be used to monitor heparin.
- Use chromogenic anti factor Xa heparin assay to monitor heparin.
 - Prophylactic range 0.1–0.4 U/mL
 - Therapeutic range 0.3–0.7 U/mL

When Do You Test for APLs?

- Low clinical suspicion
 - Venous or arterial thromboembolism \geq 50 YO
- Moderate clinical suspicion
 - Unexpected prolonged PTT without symptoms
 - Recurrent spontaneous pregnancy loss
 - Provoked venous thromboembolism (VTE) $<$ 50 YO
- High clinical suspicion
 - Unprovoked VTE or arterial thrombosis $<$ 50 YO
 - Thrombosis in unusual sites
 - Late pregnancy loss, pregnancy-related thrombosis
 - Accompanied by SLE, rheumatoid arthritis, autoimmune thrombocytopenia, autoimmune hemolytic anemia

Why detect APLs?

- Early detection of autoimmune disorder
- Explanation of venous and arterial thrombotic episodes in neurological and rheumatological conditions
- Explain recurrent spontaneous abortions, chronic low birth weights



APS Therapy

- Coumadin to prevent venous thrombosis
- Eliquis off-label (in clinical trials)
- Aspirin, Plavix, and LMWH to prevent arterial thrombosis
- Prednisone to reduce autoantibody

Bottom Line at the End (BLEAT)

The participant...

- Listed the clinical applications for a PTT mixing study.
- Listed the steps to perform a PTT mixing study.
- Explained why the mixing study is an acute care assay.
- Correlated mixing study results with lupus anticoagulant and specific inhibitor testing.



The Fritsma Factor

YOUR INTERACTIVE HEMOSTASIS RESOURCE



*Thanks for listening!
Ya got any questions?*

Sponsored by Precision BioLogic Inc
140 Eileen Stubbs Avenue
Dartmouth, NS B3B 0A9
precisionbiologic.com
Toll-free in US & CA: +1.800.267.2796