

| | APTT (S) | | |
|---|------------|-----------------|-----------------|
| Producer | IL | STAGO | SEKISUI |
| Activator | silica | silica | ellagic acid |
| Phospholipid | sythetic | rabbit cephalin | rabbit cephalin |
| Reference interval | 25.1-36.5 | 29.0-42.0 | 24.0-39.0 |
| Analyser | ACL TOP700 | STA-R EVOLUTION | CP2000 |
| Sample from patient with strong positive LA | 67.0 | 69.6 | 120.3 |
| Sample from patient with weak positive LA | 34.5 | 41.2 | 43.2 |
| sample from NOACs patient without LA | 42.9 | 51.0 | 54.7 |
| Normal pooled plasma | 31.2 | 35.3 | 31.3 |