


PrecisionBioLogic

Improving Acute Care with Coagulation Mixing Studies



George A Fritsma MS, MLS
The Fritsma Factor,
Your interactive Hemostasis Resource™
Sponsored by Precision BioLogic
Dartmouth, Nova Scotia
george@fritsmafactor.com
www.fritsmafactor.com

1

PrecisionBioLogic

Coagulation Mixing Studies Learning Objectives

The participant...

1. Defines the method and clinical applications for a PTT mixing study.
2. Lists the steps to perform a PTT mixing study.
3. Explains why the mixing study is an acute care assay.
4. Correlates mixing study results with patient coagulation testing.


2


PrecisionBioLogic

Mixing Study

An Acute Care Investigation

Purpose: *differentiate a coagulation deficiency from an inhibitor*





```
graph TD
    A[Unspecifically prolonged PT and/or APTT] --> B[Exclude pre-analytical variations  
e.g. underfilled citrate tube  
heparin contamination]
    B --> C[Perform PT and/or APTT  
mixing test]
    C --> D{ }
    D --> E[Corrections]
    D --> F[Non-corrections]
    E --> G[Perform factor assays  
if required]
    F --> H[Perform heparin and/or heparinase  
screen or specific coagulation  
factor inhibitor assays]
    G --> I[Use additional information, e.g.:  
- clinical history  
- anticoagulant therapy?  
- APTT with alternate reagent  
- Delta screen time]
    H --> I
```

Kershaw GK, Orellana D. Mixing tests: diagnostic aides in the investigation of prolonged prothrombin times and activated partial thromboplastin times. Semin Thromb Hemost 2013;39:283–90.


3

PrecisionBioLogic

Case: 32-yo Female Pre-op Screen

Six weeks post-partum

Easy bruising, frequent nosebleeds, vaginal bleeding



4

PrecisionBioLogic

Pre-op Screen

32-yo Female, 6 Weeks Post-partum

Assay	Patient	RI
HGB	11.8 g/dL	12–15 g/dL
PT	12.4 s	9.8–12.6 s
PTT (APTT)	42.5 s	25–35 s
PLT count	310,000/μL	250–450,000/μL
Fibrinogen	320 mg/dL	220–498 mg/dL
Isolated, prolonged PTT response? 1:1 PTT mix		

5

PrecisionBioLogic

Rule Out Heparin, Dabigatran

Assay	Patient	RI
TT	14 s	<21 s

- R/O unfractionated heparin (UFH) and dabigatran
 - Outpatient—consider dabigatran
 - Inpatient—unrecorded UFH flush of vascular catheter
- If dabigatran, discontinue testing, allow to clear
- If UFH, treat w/ Hepsorb (polybrene) or Hepzyme, proceed
- If no UFH, perform 1:1 PTT mix to differentiate factor deficiency from factor-specific inhibitor or “non-specific inhibitor” lupus anticoagulant (LA)

6

george@fritsmafactor.com

1


PTT Mixing Study: Cheap and Basic

- Start within 2 h to avoid specimen degradation
 - Factors V (FV) and VIII (FVIII) deteriorate
 - Platelet factors released: PF4, FV
 - Ensure patient plasma is platelet-poor, < 10,000/uL
- Mix plasma 1:1 with pooled normal plasma (NP) and perform immediate PTT on mixture
- PTT of 1:1 mix “corrects” to ≤10% above NP PTT
 - Factor deficiency
- No correction: 1:1 mix is >10% above NP PTT
 - Non-specific inhibitor, usually LA
 - Specific inhibitor (anti-FVIII), usually requires 37°C incubation

PrecisionBioLogic

7


PTT Mixing Study



PrecisionBioLogic

8

PTT Mixing Study Using 10% Rule



PrecisionBioLogic

9

1:1 PTT Mix with Incubation

- PTT of immediate mix ≤10% above NP
 - Correction: factor deficiency? But first...
 - Incubate 1:1 mix 1–2 h and repeat
- Corrected after incubated mix = factor deficiency
- No correction: PTT remains >10% above NP
 - Specific inhibitor such as anti-FVIII
 - IgG₄: Temp dependent, may require incubation
 - However, some FVIII neutralization within 10 m
 - May detect in immediate mix

PrecisionBioLogic

10

1:1 PTT Mix with Incubation


- Reflex to incubation if unincubated mix corrects
- Must also incubate normal control plasma
- Compare mix PTT to incubated normal control PTT
- May also detect temperature-dependent LA
 - ~15% of LAs are temperature-dependent

Thom J, Ivey L, Eikelboom J. Normal plasma mixing studies in the laboratory diagnosis of lupus anticoagulant. J Thromb Haemost 2003;1:2689–91

PrecisionBioLogic

11

Incubated 1:1 PTT Mix



PrecisionBioLogic

12

Mixing Study Result
32-yo Female, 6 Weeks Post-partum

Assay	Result	RI	Comment
PTT	42.5 s	25–35 s	Confirms previous PTT
PTT/control 1:1 mix immediate	32.1 s	Control 30 s	Commercial platelet-free control plasma (NP)
PTT/control 1:1 mix 1 h at 37°C	37.3 s	Control 35 s	Incubate both 1:1 mix and NP
Conclusion: both immediate <i>and</i> incubated mix PTTs correct, suspect factor deficiency, arrange for factor assays and von Willebrand disease profile			

PrecisionBioLogic


13

Factor Assay Results
32-yo Female, 6 Weeks Post-partum

Assay	Result	RI	Comment
Factor VIII	32%	50–150%	VWD?
Factor IX	92%		
Factor XI	131%		
Factor XII	113%		
HMWK	ND	65–135%	XII, HMWK & PK deficiency not associated with bleeding
PK			

PrecisionBioLogic

14



PrecisionBioLogic

15

PT and PTT Test Results
in Inherited Coagulopathies


PT	PTT	Single Factor Deficiency
Long	Normal	VII
Long	Long	X, V, II, and fibrinogen ¹
Normal	Long	VIII, IX, XI ²

¹PT & PTT prolonged when fibrinogen is <100 mg/dL, perform fibrinogen assay
²Contact factor deficiencies XII (1–3% prevalence), prekallikrein (PK, Fletcher), or high molecular weight kininogen (HMWK, Fitzgerald) also prolong PTT results, but no bleeding

PrecisionBioLogic

PTT Mix: Why Does This Work?

- Hypothetical 20% F VIII level prolongs PTT
 - PTT rgts calibrated to prolong at 30–40% FVIII, IX, XI
- Add NP with established 100% factor level
 - 1:1 mix, average of 100% and 20% = 60%, PTT corrects
- Hypothetical anti-FVIII or lupus anticoagulant
 - With typical avidity, retains ability to prolong the mix



PrecisionBioLogic

17

Case
52-yo Athletic Female

Pre-op screen for total hip replacement



18

52-yo Athletic Female
Screen Prior to Hip Replacement Surgery

Test	Result	RI
HGB	14.1 g/dL	12–15 g/dL
PT	11.2 s	9.8–12.6 s
PTT	58 s	25–35 s
PLT	170,000/ μ L	150–400,000/ μ L
Fibrinogen	410 mg/dL	220–498 mg/dL
Patient reports no bleeding or bruising, no thrombosis		

PrecisionBioLogic

19

Isolated Prolonged PTT: Differential

- Could be nothing: 5% of normals exceed limit
- Preanalytical variable: green or lavender-closure tube, hemolysis, lipemia, clotted specimen
- Outpatient: dabigatran
- Inpatient: unreported UFH
- Congenital single factor deficiency: VIII, IX, or XI, hemophilia A, B, or C with bleeding, VWD
- Congenital FXII, PK, or HMWK without bleeding
- Acquired FVIII inhibitor with severe bleeding
 - “Acquired hemophilia”
- Lupus anticoagulant (LA)



PrecisionBioLogic

20


52-yo Female PTT Mixing Study

Test	Result	Comment
TT	17 s	RI: < 21 s, rules out dabigatran
PTT	58 s	RI: 25–35 s
PTT NP	28 s	Correction if < 30.8 s (10%)
1:1 mix	35 s	25% over NP = no correction
What is the next step?		

PrecisionBioLogic

21

Acute Care Mixing Study Algorithm



PrecisionBioLogic

22

Mixing Study Considerations

- Preanalytical variables
 - Anti-Xa rivaroxaban, apixaban, edoxaban prolong PT, PTT
 - Dabigatran and UFH prolong PTT
 - Clotted, hemolyzed, lipemic specimen
 - Underfilled tube, wrong anticoagulant
 - Must be platelet-poor, <10,000/ μ L patient and NP
- Heparinase or polybrene neutralize only ≤ 1 unit/mL UFH
- Anti-FVIIIs may generate immediate neutralization
- 15% of LAs require incubation
- Weak LAs may be missed in 1:1 mix: ask for consult
 - Select a more LA-sensitive PTT reagent or prepare 4:1 mix

Mostly

PrecisionBioLogic

23

The “LA Cofactor Effect”

- Screen PTT 48 s, RI 25–35; 1:1 mix *prolongs* to 54 s
- LA “cofactor” effect thought to be prothrombin combining with LA
- LA potentiates clotting via annexin V, mix abrogates potentiation?

- Magrath M. Lupus cofactor phenomenon. Letter J Clin Pathol 1990;42:264.
- Rand JH, Wu XX, Andree HA, et al. Antiphospholipid antibodies accelerate plasma coagulation by inhibiting annexin-V binding to phospholipids: a “lupus procoagulant” phenomenon. Blood. 1998;92:1652–60.
- Clyne LP. Plasma requirement for expression of lupus-like anticoagulant. Folia Haematologica int Ma Klin Morphol Blutforsch 1986;113:841

PrecisionBioLogic

24

Normal Plasma Source?

- Home brew: ~20 normal plasmas, male ≅ female
 - Ensure plasma is platelet-poor; < 10,000/uL; PTT ≅ MRI
 - Ensure NP has ~100% of all factors
 - Elevated FVIII causes false negative results
 - Screen for LA, specific factor inhibitors. HBV, HCV, HIV
 - Aliquot and freeze
- Or purchase commercial plasma
 - GMP & frozen meets all criteria
 - Lyophilized acceptable when validated
 - Processed with stabilizers



Clinical and Laboratory Standards Institute. One-stage prothrombin time (PT) test and activated partial thromboplastin time test (APTT) approved guideline—second edition. CLSI Document H47-A2. CLSI, Wayne PA. 2008.

25

What Limit Defines Correction?

No Consensus; Fritsma Factor Quick Question Results

- Limits based on fixed PTT value
 - 1:1 mix within RI upper limit (95% or 99% CI, **39%**)
 - 1:1 mix within RI upper limit + 5 seconds (**8%**)
- Limits based on NP PTT value
 - 1:1 mix within NP PTT value + 5 seconds (**14%**)
 - 1:1 mix within NP PTT + 10% (same as ratio 1:1, **32%**)
- Limit formula using patient, NP, and 1:1 mix
 - **Must incubate patient sample, NP, and 1:1 mix**
 - Chang's % deviation; Rosner index
- Other (**7%**): combination of RI and Rosner
 - Dedicated RI for mix

PrecisionBioLogic

26

1:1 Mix Limit Based on % Correction

$$\% \text{ Correction} = \frac{\text{Patient PTT} - 1:1 \text{ mix PTT}}{\text{Patient PTT} - \text{NP PTT}} \times 100$$

$$\% \text{ Correction} = \frac{42.5 - 32.1 = 10.4}{42.5 - 30 = 12.5} = 0.83 = 83\%$$

$$\text{Factor Deficiency} = \geq 75\%$$

$$\text{Inhibitor} = < 75\%$$

% Correction verified by local laboratory

Chang SH, Tillema V, Scherr D. A "percent correction" formula for evaluation of mixing studies. Am J Clin Pathol 2002;117:62-73.

PrecisionBioLogic

27

Rosner Index

$$\text{Rosner Index} = \frac{1:1 \text{ mix PTT} - \text{NP PTT}}{\text{Patient PTT}} \times 100$$

$$\text{Rosner Index} = \frac{32.1 - 30}{42.5} \times 100 = 4.9$$

$$\text{Inhibitor} \geq 11$$

$$\text{Correction} < 11$$

Rosner index validated by local laboratory

Rosner E, Pautzner R, Lusky A, Modan M, Many A. Detection and quantitative evaluation of lupus circulating anticoagulant activity. Thromb Haemost 1987; 57: 144-147.

PrecisionBioLogic

28

PrecisionBioLogic

59-yo Male Former Hockey Player

Total knee replacement preop



29

59-yo Male Former Hockey Player Screen Prior to Knee Replacement Surgery

Test	Result	RI
HGB	14.8 g/dL	12–15 g/dL
PT	11.2 s	9.8–12.6 s
PTT	38 s	25–35 s
PLT	310,000/ μ L	150–400,000/ μ L
Fibrinogen	390 mg/dL	220–498 mg/dL
Patient reports no bleeding or bruising, no thrombosis		

PrecisionBioLogic

30

When to Perform Mixing Study

- Any PTT > RI upper limit
- Any PTT > RI upper limit + 5 seconds
- Any PTT > RI upper limit with consult
 - Is patient bleeding or clotting?
 - Possible “weak” LA: use 4:1 mix
 - Lupus sensitive PTT reagent
 - Factor sensitive PTT reagent

Pengo V, Tripodi A, Reber F, et al. Update of the guidelines for lupus anticoagulant detection. J Thrombos Haemost 2009;7:1737–40.

PrecisionBioLogic

31

59-yo Male Former Hockey Player

Test	Result	Comment
TT	17 s	RI: < 21 s, rules out dabigatran
PTT	38 s	RI: 25–35 s
PTT NP	31 s	Correction if < 34.1 s (10%)
1:1 mix	35 s	Correction? No correction?
What is the next step?		

PrecisionBioLogic

32

59-yo Male Former Hockey Player Clinical Consult

- Consult: if he is well, go on to procedure
- Prior thrombotic events (VTE)
 - Perform mix using 4:1 patient plasma to NP
 - Or choose PTT reagent that is LA-sensitive
- If anatomic bleeding symptoms, test FVIII, FIX, FXI
 - Vitamin K deficiency
 - Renal insufficiency
 - Liver disease, malignancy, VWD

PrecisionBioLogic

33

PrecisionBioLogic

2-yo Hemophilic Boy



34

2-yo Hemophilic Boy

Test	Result	RI
HGB	11.8 g/dL	9.6–15.6 g/dL
PT	11.2 s	9.8–12.6 s
PTT	65 s	25–35 s
PLT	310,000/ μ L	150–400,000/ μ L
Fibrinogen	390 mg/dL	220–498 mg/dL
Inflamed, swollen knee and ankle		

PrecisionBioLogic

35

Mixing Study Result 2-yo Hemophilic Boy


Assay	Result	RI	Comment
PTT	65 s	25–35 s	Confirms previous PTT
PTT/control 1:1 mix immediate	33.5 s	Control 30 s	
PTT/control 1:1 mix 1 h at 37°C	47.9 s	Control 35 s	Control is incubated alone and with mix
Conclusion: Anti-FVIII inhibitor			

PrecisionBioLogic

36

Factor VIII Assay

- Dilute plasma 1:10
- Add factor VIII-depleted reagent plasma 1:1
- Add PTT reagent, incubate 3 minutes
- Add CaCl₂, record interval to clot formation
- Compare result in seconds to calibration curve



Factor VIII Assay Dilutions Parallelism Indicates No Inhibitor

Plasma Dilution	Seconds	Raw Factor VIII Activity	Computed Factor VIII Activity (x dilution)
1:10 "undiluted"	90 s	20%	20%
1:20	104 s	10%	20% (parallel)*
1:40	107 s	5%	20% (parallel)
1:80	110 s	2.5%	20% (parallel)

* <10% difference from undiluted indicates parallelism, no inhibitor

FVIII Assay Dilutions non-Parallelism Indicates Inhibitor

Plasma Dilution	Seconds	Raw Factor VIII Activity	Computed Factor VIII Activity (x dilution)*
1:10 "undiluted"	80 s	10%	10%
1:20	93 s	8%	16%
1:40	107 s	5%	20%
1:80	108 s	4%	32%

* >10% difference from undiluted, rising = non-parallel, implies inhibitor

Kasper CK. Laboratory diagnosis of factor VIII inhibitors. In Kessler C, Garvey MB, Green D, Kasper C, Lusher J. Acquired Hemophilia 2nd Edition. Excerpta Medica 1995

55-yo Male with Atrial Fibrillation

Test	Result	RI
HGB	13.8 g/dL	12–15 g/dL
PT	17.2 s	9.8–12.6 s
PTT	159 s	25–35 s
PLT	310,000/ μ L	150–400,000/ μ L
Fibrinogen	20 mg/dL	220–498 mg/dL

55-yo Male with Atrial Fibrillation

Assay	Result	RI
PTT	159 s	25–35 s
TT	> 150 s	< 21 s
PTT/control 1:1 mix immediate	78 s	Control 30 s
PT/control 1:1 mix immediate	15.2 s	Control 12 s
What do you recommend?		

If the PT is Prolonged

- Congenital deficiencies of II, V, VII, or X
 - PT and PTT long: II, V, X
 - PT only: VII, skip mixing and go to factor assay
 - Prevalence: 500,000–1:2,000,000
- Liver disease: PT prolongs before PTT due to des-carboxy II, VII, and X, reduced factor V
- Vitamin K deficiency: des-carboxy II, VII, and X
- Anti-Xa direct oral anticoagulants
 - Rivaroxaban, apixaban, edoxaban

Isolated Prolonged PTT: Summary

- Random benign prolongation, 95% CI
- Lupus anticoagulant: Prevalence of 1–3%
- Drug reaction producing transient LA
- Unrecorded heparin, dabigatran, oral anti-Xa
- Known hemophilic who fails FVIII concentrate Rx
- Hemorrhage or ecchymoses signal acquired coagulopathy; vitamin K deficiency, liver disease
- Specific inhibitor, anti-FVIII
 - Postpartum, malignancy
 - Autoimmune disorders, > 60-yo

Sahud MA. Factor VIII inhibitors. Laboratory diagnosis of inhibitors Semin Thromb Hemost 2000;26:195–203.
PrecisionBioLogic

43

Develop Mixing Study Reliability

- PTT reagent sensitivities
 - 30–40% FVIII, FIX, FXI
 - Intermediate sensitivity to LA
- NP consistency: ~100% activity for all factors
- Consultation for equivocal patient results
- Employ consistent correction limit

PrecisionBioLogic

44

Perform Mixing Studies Locally

- Unexpected isolated prolonged PTT or PT may require immediate action
- Local results may immediately direct therapy
- Delayed specimen may deteriorate
- Forward results to ref lab to direct follow-up

PrecisionBioLogic

45

PrecisionBioLogic

BLATE
(Bottom Line at the End)

1. Define the method and clinical applications for a mixing study.
2. List the steps to perform a mixing study
3. Explain why the mixing study is a first-line assay.
4. Correlate mixing study results with patient coagulation testing.

Precision BioLogic Inc
140 Eileen Stubbs Avenue
Dartmouth, NS B3B 0A9
precisionbiologic.com
Toll-free in US & CA: +1.800.267.2796

46

PrecisionBioLogic

Thanks for
listening!



Questions?

The Fritsma Factor,
Your interactive Hemostasis Resource
Sponsored by Precision BioLogic
Dartmouth, Nova Scotia
george@fritsmafactor.com
www.fritsmafactor.com

47